



## Simplicity

NO calibration required, easy to use.

## Portability

Lightweight, handheld.

## Performance

Specifically calibrated and validated for dogs and cats.

**aLcose**<sup>®</sup>

## Specification

- Enzyme: Glucose Dehydrogenase (GDH-FAD)
- Sample Type: Whole Blood (Capillary or Venous)
- Measurement Time: 4 Seconds
- Measurement Range: 20-600 mg/dl (1.1~33.3 mmol/l)
- Blood sample volume: 0,4ul
- Date Storage: 400 Result with date and time
- Percision: CV < 3%
- Hematocrit Range: 10-70%



1 Insert strip and auto turn on meter.




2 Apply a drop of blood sample.



3 Obtain result in 4 seconds.

## Summary of Accuracy Results Relative to IDEXX

● Cats



● Dogs

